

HEM-618 HbA1c Meter

Technical specifications

- ▶ Test parameters: HbA1c, eAG
- ▶ Sample type: whole blood
- ▶ Sample volume: 10µL
- ▶ Test time: 3-4 min.
- ▶ Measurement range: HbA1c 4.0%-15.0%
- ▶ CV: ≤5%
- ▶ Storage: 10,000 sample; 12 months QC value, power failure protection
- ▶ COM: RS-232, HIS/LIS
- ▶ Power requirement: AC220V±10%, 50Hz,50VA, with Li battery option
- ▶ Dimension: 205mm×130mm×70mm, 1KG
- ▶ Working temperature: 5-40°C
- ▶ Consumables: HbA1c cartridge, kits, IC card



Auf der Heide 15 - D-37351 Dingelstädt
Federal Republic of Germany

A + A Diagnostics

CE
ISO13485

HEM-618 HbA1c Meter



- Dry system, maintenance free
- Multi-color touch screen
- Disposable cartridge, point of care